



Single-use RhinoLaryngoscope

RH Series



Single-use



Portable



No Cross-infection



Always Ready To Go

—

For your easier
Rhinolaryngoscopy
in OPD



MD722076

An Endoscopy Solution Company


P/N: ENG-EU-CATA-SPE-fRHS-suRH-INTENT-202310V1.0

Specification

Model	O.D.	I.D.	Deflection	Field of View	Working Length	Direction of View
RH-M22	2.2mm	0mm	210° Up/ Down	120°	350mm	0°
RH-M32	3.1mm	1.2mm	210° Up/ Down	120°	350mm	0°
RH-M40	4.0mm	2.0mm	210° Up/ Down	120°	350mm	0°
RH-M50	5.0mm	2.8mm	210° Up/ Down	120°	350mm	0°
RH-M52	5.0mm	2.6mm	210° Up/ Down	120°	350mm	0°
RH-M58	5.8mm	3.4mm	210° Up/ Down	120°	350mm	0°

RH Series Imaging Part

Image processing, display and recording
All-in-one



Specification	
Model	MS-8
Dimension	386.4 (h) x 253.6 (w) x 38.5 (d) mm
Weight	2.6 kg
Display size	15.6 inches
Display resolution	1080P(1920*1080)px
Display type	TFT-LCD, touch screen
Function	Dual view, WIFI, Brightness control, Contrast control, Using time calculator, etc.
Working time	Continuous working time ≥4 h
Memory capacity	Built-in SD card standard 64 GB

Ordering information

Product	Item code
Single-use RhinoLaryngoscope	
RH-M22	2018B-PA005
RH-M32	2018B-PA004
RH-M40	2018B-PA003
RH-M50	2018B-PA001
RH-M52	2018B-PA002
RH-M58	2018B-PA006
Image processor	
MS-8	210-0697-00

Specifications, design and accessories are subject to change without any notice or obligation on the part of the manufacturer.

HugeMed

Shenzhen HugeMed Medical Technical Development Co., Ltd.
415, 416-1, 516-1, Building 2, No.1, Mawu Road, Baoan Community, Yuanshan Street,
Longgang District, Shenzhen, Guangdong, 518115, China



+86-755-22275833



+86-755-22275833



info@hugemed.net

For more information
please visit our website:
www.hugemed.net